Candulor AG

CH-8602 Wangen / ZH, Pünten 4
Tel. +41 (0)44 805 90 00, Fax +41 (0)44 805 90 90
candulor@candulor.ch, www.candulor.com

Subsidiary: Candulor USA Inc.

7462 N. Figueroa St., Suite 104, Los Angeles CA 90041
Phone +1 (323) 254-1430, Phone +1 800 436-3827
Fax +1 (323) 254-5146
info@candulor.us, www.candulor.com

Candulor Composite

Mikrogefülltes Verbundmaterial
Micro filled bonding agent
Matériau de liaison microchargé
Materiale de unione microriempito
Bonding micro-relleno
Microgevuld tandheelkundig hechtmiddel

Verarbeitungsanleitung
Instructions for Use
Instructions de mise en oeuvre
Istruzioni d'uso
Instrucciones de uso
Gebruiksaanwijzing

Candolor Composite

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Instructions for Use

Description

Candolor Composite is a micro-filled bonding agent designed to achieve an optimum bond between two different dental resins or composites.

Indication

Candolor Composite is used as bonding material for mould modifications and colour characterisations of Candolor teeth made of PMMA, TCR or NFC.

Colour characterisations are performed with:

- Candolor Stains for Resin Teeth

Mould modifications are performed with:

- SR Adoro Ivoclar Vivadent
- SR Chromasit Ivoclar Vivadent

Contraindications

- Complementing of ceramic-based denture teeth
- Do not use Candolor Composite if the patient is known to be allergic to any of the product's ingredients.

Side effects

In rare cases, an allergic reaction may occur. Candolor Composite should not be used in patients with suspected or known allergy to any of its ingredients.

Warning

Candolor Composite is designed for use in dental technology only. Avoid contact of uncured material with the eyes and the skin. Contact with uncured Candolor Composite may cause slight irritation of the skin and induce sensitization to methacrylates. Many of the gloves commonly used do not provide effective protection against the sensitizing effect of methacrylates.

Interaction

None known to date.

Storage

- Store packages at room temperature (12–28 °C/54–82 °F)
- Close syringes immediately after use (exposure to light may lead to premature polymerization)
- Do not use the material after the shelf life expiry date
- Keep out of the reach of children

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

Application

1. Pre-treating the contact surface

Roughen the resin surface with a cross-cut carbide cutter to achieve an optimum bond. Produce the desired preparations, such as grooves, for characterization. Blast tooth surface with abundant aluminium oxide, grain size 50–100 µm pressure. Keep the degreased surface clean. Thoroughly wet the clean surface with denture base monomer (e.g. Candolor Aesthetic). Allow the material to react for approx. 2–4 minutes.

Important

It is mandatory both to wet the contact areas thoroughly with denture base monomer and allow the monomer to react in order to avoid weakening of the bond. Keep surface clean.

2. Application of Candolor Composite

Extrude Candolor Composite from the syringe directly onto the contact surface. Evenly spread the material on the surface with a disposable brush or a spatula.

Important

Observe layer thickness of 0.2 to 0.5 mm.

3. Working time

The material takes approx. 3 minutes to cure when exposed to a luminous intensity of 3000 lux.

4. Polymerization

Polymerization may be carried out under direct light with the following units:

Candolor recommends the following light curing units

Apparatus	Manufacturer	Recommended curing time
Quick	Ivoclar Vivadent	60 seconds
Lumamat 100	Ivoclar Vivadent	11 minutes (P2)
Spectramat	Ivoclar Vivadent	4 minutes

After polymerization, the surface will be slightly sticky and must therefore be protected from dirt. The sticky layer facilitates the subsequent processing of veneering materials.

5. Finishing

After polymerisation, the teeth are finished in accordance with the instructions for use for Stains for Resin Teeth (Candolor) or for SR Adoro (Ivoclar Vivadent) and SR Chromasit (Ivoclar Vivadent). If finishing is not carried out immediately after polymerization, store reconstruction in a clean, dark place.

Important

Check shadow areas for thorough polymerization.

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